


# VISICOIL™

LINEAR FIDUCIAL MARKER

new!



## Designed for Interventional Radiology Approved by Interventional Radiology

VISICOIL's 0.50mm and 0.35mm markers are now available pre-loaded in 21 and 22 gauge 8cm and 15cm needles. Ideal for Interventional Radiology allowing for easier access with improved patient comfort and safety due to the smaller needle size.

### IMPLANTATION SAFETY

- Reduced risk of pneumothorax
- Reduced risk of bleeding / edema
- Improved safety and comfort for patient

### IMPLANTATION EFFICIENCY

- Shorter 8cm and 15cm needles ideal for standard bore CT
- Steerable needle aids in accurate marker placement
- Centimeter marks for depth placement



[www.visicoil.com](http://www.visicoil.com)

Protect,  
enhance  
and save  
lives

Iba

# Safety, Efficiency and Clinical Confidence in Tumor Targeting

### PATIENT SAFETY

- Conventional gold markers require a 17g needle
- VISICOIL IR needles are up to 50% smaller




- 17g
- 21g
- 22g

VISICOIL Needle Size

### TREATMENT EFFICIENCY

- One marker - multiple reference points
- Flexible marker - visualize tissue deformation & motion



Pancreas: 2 VISICOIL 0.50mm x 1cm

### CLINICAL CONFIDENCE

- Minimized artifact for accurate tumor targeting
- Non-migrating design



Liver: 1 VISICOIL 0.35mm x 1cm

### NEW VISICOIL IR PRODUCT RANGE

VISICOIL SIZE	NEEDLE SIZE	PART NUMBER
0.50mm x 0.5cm	21 gauge x 8cm	SN-050-005-PL08
0.50mm x 0.5cm	21 gauge x 15cm	SN-050-005-PL15
0.50mm x 1.0cm	21 gauge x 8cm	SN-050-005-PL08
0.50mm x 1.0cm	21 gauge x 15cm	SN-050-010-PL15
0.35mm x 0.5cm	22 gauge x 8cm	SN-035-005-PL08
0.35mm x 0.5cm	22 gauge x 15cm	SN-035-005-PL15
0.35mm x 1.0cm	22 gauge x 8cm	SN-035-010-PL08
0.35mm x 1.0cm	22 gauge x 15cm	SN-035-010-PL15

VISICOIL is also available loose in various sizes and lengths to be used in the needle of your choice.

### CONTACT DETAILS

#### NORTH AMERICA | LATIN AMERICA

RadioMed Corporation  
3150 Stage Post Drive, Suite 110  
Bartlett, TN 38133  
Tel: +1 901 432 7202  
Fax: +1 901 432 7206  
Toll Free: 866 429 0922  
visicoil@iba-group.com

#### EUROPE | MIDDLE EAST | ASIA PACIFIC

IBA Dosimetry GmbH  
Bahnhofstr.5  
90592 Schwarzenbruck, Germany  
Tel: +49 9128 607 0  
Fax: +49 9128 607 10  
email: info@iba-dosimetry.com

### PLEASE SEE BELOW CODE FOR IMPLANTATION INSTRUCTIONS



CALL CUSTOMER SERVICE FOR YOUR LOCAL REPRESENTATIVE  
800 429 0922